

-86°C Ultra Low Temperature Freezers



ULTRA LOW TEMPERATURE FREEZERS

The Ideal – 86°C Freezing Environment in Capacities from 17 cu.ft. to 24 cu.ft.

Ideally suited for ultra storage laboratories, for long-term preservation and storage of blood, specimens and components, and in testing of various types. Ranging in size from economical chest freezer to large capacity upright configuration, one of these models is sure to be suitable for your needs. Whatever your preservations needs are, SANYO provides the right equipment to meet your exact needs. SANYO preservation systems employ advanced technology to insure a high precision temperature environment.

> innovation performance reliability support



preservation

-86°C PRODUCT FAMILY

A wide range of Ultra-low temperature freezers.

This product family offers different types and models of ultra-low freezer which can achieve and maintain temperatures down to -86°C (-123°F) at ambient temperature of 30°C (85°F). It is a wide selection, offering a variety of sizes and types suitable for uses such as long term preservation, testing, frozen storage etc.

Specially Designed Compressor for Ultra-Low Temperature.

The refrigeration system is at the heart of producing ultra-low temperature freezing. In the process of developing various technical solutions in order to achieve and maintain ultra-low temperatures, Sanyo designed a compressor especially for ultra low-temperature use. Sanyo is a leading manufacturer of ultra-low temperature freezers that develops and produces all its own components, from the refrigeration units to the microprocessors.


An Ideal Freezing Environment By Means of A Double-Insulation Wall and the "Hot Line".

In Ultra-Low temperature freezers, effective insulation is extremely critical since the temperature difference between the outside and the inner compartment



DESIGN

can be as high as 115°C. The insulation used in Sanyo's ultra-low temperature freezers is two layers of rigid, foamed in place polyurethane. This design prevents the insulation layers from distortion and cracking that might occur due to temperature differences inside and out, and creates the most efficient insulation material available today. Moisture condensation at the top edges of the cabinet due to differences in temperature inside and out causes frost and icing problems, which may reduce heat insulation efficiency and obstruct door movements. They are prevented by the "hot line" by means of which hot gas from the higher temperature circuit is circulated through the problem areas.






UPRIGHT -86°C ACCESSIBILITY

Evolutionary Design

• The newly developed refrigeration system and freezer structure minimizes intrusive operating noise.

The newly developed back cover combined with new aerodynamically designed and placed components in the refrigeration compartment provide superior air flow, making it possible to drastically reduce the stress to the freezer and contributing to excellent durability.

Two independent and insulated inner doors ensure maximum interior chamber uniformity.



Safety

■ High and low temperature warning provides an audible and visual alarm when the temperature deviates more than ±5°C to ±20°C (adjustable) from the set point.

Alarm ring back function ensures that buzzer will resume operation, should alarm conditions continue after it is silenced.

 Microprocessor-controlled filter-clog check function

protects the refrigeration circuit.

• Control panel with digital display for easy operation.





preservation

Accessibility

• A new rugged inner door latch tightly closes the inner door and makes accessing chamber contents easier.

• The condenser filter is located at the bottom right side of the front panel to make filter removal and cleaning easier.

• Alarm lamp and buzzer offer secure warning of power failure or abnormal temperature increase.

New rugged, one-handed outer door latch has a hasp to allow a padlock for additional security





SAFETY & ALARM

Built-In Temperature & Power Failure Alarms(Lamp/Buzzer).

In case of power failure or an irregular rise in temperature, a rechargeable, Cadnica battery-operated indicator lamp and alarm will be activated. Optional CO_2 and LN_2 back-up systems are available and are self-activated when a power outage occurs. This equipment helps ensure that the contents will be protected in the event of any power failure or mechanical trouble.

Easy to Operate.

Overall operability and dependability are greatly enhanced by improved design details. Highly durable hinges are used to support the large door so both opening and closing are facilitated. Door handles are equipped with a latch-locking system.

A Wide Range to Choose From Large, Small, Chest and Upright- A freezer to suit your needs.

Chest type freezers low profile design makes the placement and removal of contents very convenient. The upright type is designed to accommodate a large capacity in a limited space. Even very narrow spaces can be readily utilized.

RELIABILITY & CONTROL

A Special Refrigerant and Oil to Maintain Stability and Reliability.

In order to expand refrigeration capacity, a special compound refrigerant has been used. This refrigerant relieves the refrigeration system of much of its load by evaporating at a lower temperature level within the circuit. Also, a special grade of refrigeration machine oil is used that offer superior properties of high rate of recovery and outstanding resistance to heat and wear. This contributes to the long life and dependability of the freezer unit as a whole.

Microprocessor Temperature Control With Digital Design.

The temperature inside the freezer can be set and monitored easily by means of precise microprocessor temperature control with digital display. The thermostat utilizes a platinum resistor (Pt 100) sensor which is precise and extremely durable.

MDF-593C



17.1 cu.ft.

MDF-793C



Specifications

Madal	MDF-593C	MDF-793C	
Model			
Temperature range	-20°C to -86°C	-20°C to -86°C	
Exterior dimensions	2,010 x 770 x 1,070	2,570 x 770 x 1,070	
(W x D x H)	(79.1 x 30.3 x 42.1)	(101.2 x 30.3 x 42.1)	
Interior dimensions	1,280 x 500 x 762	1,840 x 500 x 762	
(W x D x H)	(50.3 x 19.7 x 30)	(72.4 x 19.7 x 30)	
Effective capacity	487 liters (17.1 cu.ft.)	701 liters (24.7 cu.ft.)	
Exterior Cabinet	Painted Steel		
Interior Cabinet	Stainless Steel		
Inner Lid	3	4	
Insulation	Rigid polyurethane foamed-in place		
Compressor High stage si	de Hermetic type 750W	Hermetic type 1,100W	
Low stage si	e Hermetic type 750W	Hermetic type 1,100W	
Evaporator High stage si	Cascade condenser		
Low stage si	Tube on sheet (shared with interior)		
Condenser High stage si	Fine and tube type		
Low stage si	e Shell and tube	Shell and tube type	
	Microprocessor: Keypad input		
Temperature Control	Temp. input range: -20°C to -95°C (1C increment)		
	Set value memory: non-	volatile memory	
Temperature Display	Digital Display		
Sensor	Platinum resistanc	Platinum resistance (Pt. 100)	
Alarm System	Selectable high temp. alarm (+10°C & +15°C fro	Selectable high temp. alarm (+10°C & +15°C from set point), Power Failure alarm, Filter	
	check lamp, Remote alarm contact		
Net Weight (Approx.)	305kg (672 lbs)	375kg (827 lbs)	
Power Requirement	220V / 60hz Built-in voltage booster		

*Appearance and specifications are subject to change without notice.

Model	MDF-U7386SC	MDF-U5386SC
Temperature range	-50°C to -86°C (1°C increment)	-50°C to -86°C (1°C increment)
Maximum cooling	-86°C	-86°C
performance	(Ambient temp. 30°C)	(Ambient temp. 30°C)
Exterior dimensions	44.5" x 34.4" x 78.3"	35.0" x 34.4" x 78.3"
(W x D x H)	(1130 x 875 x 1990mm)	(890 x 875 x 1990mm)
Interior dimensions	34.2" x 23.6" x 50.4"	24.8" x 23.6" x 50.4"
(W x D x H)	(870 x 600 x 1280mm)	(630 x 600 x 1280mm)
Net weight	783 lbs. (Approx. 355kg)	672 lbs. (Approx. 305kg)
Effective capacity	23.5 cu.ft. (668 L)	17.1 cu.ft. (483 L)
Shelf	Stainless steel, Adjustable,	Stainless steel, Adjustable,
	3 shelves, W848 x D533mm,	3 shelves, W608 x D533", Max.
	Max. load: 50kg (110 lbs.)/shelf	load: 50kg (110 lbs.)/shelf
Access port	17mm diameter, 3 locations	.07"(17mm) diameter, 3 locations
	(back, bottom left/right corner)	(back, bottom left/right corner)
Compressor	Hermetic type, Output: 1100 W	Hermetic type, Output: 1100 W
	(high stage), 1100 W (low stage)	(high stage side), 1100 W (low
		stage side)
Refrigerant	HFC refrigerants	HFC refrigerants
Alarm	High/low temperature,	High/low temperature,
	Power failure, Filter check,	Power failure, Filter check,
	Battery check	Battery check
Remote alarm contact	Allowable contact capacity:	Allowable contact capacity:
	DC 30V, 2A	DC 30V, 2A
Accessories	1 set of keys, 1 scraper	1 set of keys, 1 scraper
Power Requirement	220V / 60 Hz	220V / 60 Hz/ 1 phase
	Built-in voltage booster	Built-in voltage booster
		Nema 6-15 Recepticle



Rack Capacity: 24





Strip Chart Recorder MTR-85H

> Recorder mounting bracket **MDF-S3085** (for MTR-85H)

CVK-UB2 (I) LN2 backup system CVK-UBN2

ISO 9001 & 14001 certified



SANYO Commercial Solutions A Division of SANYO North America 1062 Thorndale Avenue, Bensenvile, L60106 USA Toll Free USA 800-858-8442 • Fax 630-238-0074 www.sanyobiomedical.com SANYO Canada, Inc. 1- 300 Applewood Crescent, Concord, Ontario L4K 5C7 905-760-4025 • Fax 905-760-9405

SANYO

www.sanyobiomedical.com

σ